



PROTEAM FLEXFIT® BACKPACK VACUUM HARNESS PROMOTION

May 1, 2016 – June 30, 2016

Purchase any FlexFit articulating backpack harness to upgrade existing applicable vacuum models and receive the below SPECIAL discounts.

FlexFit Backpack Harness Promotion Discounts

# of Units	% Discount	Promo Code
1 - 4	49%	PT162FFH1
5 +	51%	PT162FFH2

Qualifying Equipment:

840011 – FlexFit Backpack Harness Assembly

The FlexFit harness fits all of the below ProTeam backpack vacuum models:

GoFree® Flex Pro battery

Super Coach Pro® 6 & 10

ProVac® FS 6

- Only the listed items qualify for the above promotional discount.
- Quantity discounts apply to single orders shipped to one location.
- **PROMO CODE(S): MUST BE NOTED ON ORDER.**
- Discounts are off of list price.
- Standard account freight terms apply.
- Distributor may include other ProTeam vacuums, filters & accessories to an order, as applicable, to help qualify for free ground freight based on their account terms.
- *Bids, Contract Pricing and Deviated Programs are excluded.*

***All orders must be received and scheduled for immediate shipment between
May 1, 2016 – June 30, 2016, to receive the promotional discount.***

Please contact ProTeam Customer Service toll-free at 1-866-888-2168, by fax at 800-844-4995, or via email at customerservice.proteam@emerson.com, or contact your local ProTeam Representative, to place your order.

This message was sent as part of our existing business relationship. To update your contact info or if you'd prefer to be removed from further distribution lists, please email customerservice.proteam@emerson.com with the correct info or with your removal request.