



INTERNAL USE ONLY

CONFIDENTIAL

For ProTeam and ProTeam Manufacturer Rep Information Only

Representative Demo Equipment Pricing

Updated 01/24/17

| Machine Model | SKU # | List price | Under 1 year, slightly used condition | 1 to 2 years old, normal demo use and condition | 1 year or older with HEAVY USAGE, or 2 or more years old |
|---------------|-------|------------|---------------------------------------|---|--|
|---------------|-------|------------|---------------------------------------|---|--|

BACKPACKS, CANISTERS AND HIP-STYLE VACS

| | | | | | |
|----------------------------|--------|----------|----------|----------|----------|
| AviationVac | Any | \$619.00 | \$247.60 | \$216.65 | \$185.70 |
| MegaVac | Any | \$559.00 | \$223.60 | \$195.65 | \$167.70 |
| ProVac FS 6 (100727 kit) | 107363 | \$549.00 | \$219.60 | \$192.15 | \$164.70 |
| ProVac FS 6 (Wessel kit) | 107461 | \$889.00 | \$355.60 | \$311.15 | \$266.70 |
| ProVac CN canister | Any | \$629.00 | \$251.60 | \$220.15 | \$188.70 |
| RunningVac canister | Any | \$589.00 | \$235.60 | \$206.15 | \$176.70 |
| LineVac ^{er} HEPA | Any | \$749.00 | \$299.60 | \$262.15 | \$224.70 |
| QuietPro BP HEPA | Any | \$609.00 | \$243.60 | \$213.15 | \$182.70 |
| QuietPro CN HEPA canister | Any | \$659.00 | \$263.60 | \$230.65 | \$197.70 |
| Super HalfVac Pro | Any | \$529.00 | \$211.60 | \$185.15 | \$158.70 |
| Super CoachVac HEPA | Any | \$539.00 | \$215.60 | \$188.65 | \$161.70 |
| Super QuarterVac HEPA | Any | \$539.00 | \$215.60 | \$188.65 | \$161.70 |
| Super Coach Pro 10 | Any | \$539.00 | \$215.60 | \$188.65 | \$161.70 |
| Super Coach Pro 6 | Any | \$529.00 | \$211.60 | \$185.15 | \$158.70 |

UPRIGHT MODELS

| | | | | | |
|----------------------|--------|----------|----------|----------|----------|
| ProForce 1200XP HEPA | 107251 | \$459.00 | \$183.60 | \$160.65 | \$137.70 |
| ProForce 1500XP HEPA | 107252 | \$569.00 | \$227.60 | \$199.15 | \$170.70 |
| ProGen 12 | 107329 | \$489.00 | \$195.60 | \$171.15 | \$146.70 |
| ProGen 15 | 107330 | \$539.00 | \$215.60 | \$188.65 | \$161.70 |

BATTERY MODELS

| | | | | | |
|-----------------------|-----|------------|----------|----------|----------|
| GoFree Flex Pro (NEW) | Any | \$1,699.00 | \$679.60 | \$594.65 | \$509.70 |
|-----------------------|-----|------------|----------|----------|----------|

WET/DRY MODELS

| | | | | | |
|-------------------------|--------|----------|----------|----------|----------|
| ProGuard LI 3 Cordless | 107495 | \$339.00 | \$135.60 | \$118.65 | \$101.70 |
| ProGuard 4 Portable | 107128 | \$239.00 | \$95.60 | \$83.65 | \$71.70 |
| ProGuard 10 | 107129 | \$509.00 | \$203.60 | \$178.15 | \$152.70 |
| ProGuard 15 | 107130 | \$759.00 | \$303.60 | \$265.65 | \$227.70 |
| ProGuard 15 w/ squeegee | 107359 | \$949.00 | \$379.60 | \$332.15 | \$284.70 |
| ProGuard 16 MD | 107386 | \$499.00 | \$199.60 | \$174.65 | \$149.70 |
| ProGuard 20 | 107131 | \$809.00 | \$323.60 | \$283.15 | \$242.70 |
| ProGuard 20 w/ squeegee | 107360 | \$999.00 | \$399.60 | \$349.65 | \$299.70 |
| ProBlitz | 107132 | \$289.00 | \$115.60 | \$101.15 | \$86.70 |
| ProBlitz XP | 107133 | \$329.00 | \$131.60 | \$115.15 | \$98.70 |

WARRANTY COVERAGE NOTE: Limited Warranty coverage on vacuums is restricted to the original owner/purchaser. In the case of Manufacturer Rep demos, however, ProTeam will make an exception and any time remaining from the date of original invoice can be extended to a buyer of a Rep Demo machine. Example, a 2-year old machine with 1 year remaining on the Limited Warranty can have that remaining 1 year of coverage included in the sale price. Note that no warranty coverage exists on accessories.